

ManaWaterSling

Waterproof Arm Sling

Part No: MWSLING01



INSTRUCTIONS FOR USE

Before using, please read the following instructions completely and carefully. Correct application is important to the proper functioning of this sling.

Indications for Use

The ManaWater Sling is designed to be worn in the shower, bath, pool or beach to help in the healing process for a variety of shoulder injuries and surgeries including:

- Dislocated Shoulder
- Rotator Cuff Tear and/or Repair
- Clavicle or Humerus Fracture
- Neuropathy and other Nerve Issues
- Fractured Arm
- AC Separation
- Tendonitis

Usage Precautions

Always check the condition of the sling before use. Application should be done with great care to ensure efficiency. For hygiene reasons, it is proposed to be used by a single patient. In exceptional cases of very sensitive skin, irritation may be caused. In this case, please contact your physician.

Application:

1. Position wrist loop in front of body
2. With uninjured arm and shoulder, slide large cord strap over head.
3. Position cushioned neck/shoulder strap where comfortable.
4. Place hand of injured shoulder into the wrist loop.

User Maintenance, Storage, and Cleaning

1. Check the condition of the sling before use.
2. Hand wash with warm or cold water. Use a soft detergent. Do not use bleach.
3. Air dry Do not wash in a washing machine or dry in a tumble dryer.



ManaMed Inc.
5240 W Charleston Blvd.
Las Vegas NV 89146
United States
Ph: 888-508-0712
www.ManaMed.com



MANAMED™

©2021 ManaMed Inc. All rights reserved.
ManaMed, the ManaMed logo,
ManaWaterSling logo,
are trademarks of ManaMed Inc.
Registered in the U.S.