

### Features:

Excellent choice for patients who are unable to pull a knee brace over their feet.

Medial and lateral hinge straps provide excellent support.

Unique wrap design provides for an excellent custom fit.

### How to Apply:

Unfasten all straps and closures. Position the back with the popliteal opening and behind the knee so popliteal opening falls at the crease behind the knee. Align the hinges with the center at the knee joint. Wrap around the front and secure the hook and loop closure. Support should fit snug but not so tight as to restrict circulation. Hinges should be aligned with the prominent bones on each side of the knee. Secure straps beginning with top strap and working your way downward. Feed the straps through the rings and secure in place.

### Setting the Hinge:

Hinge has optional flexion and extension settings. The caregiver should determine the correct setting for the condition. Remove hinges from pockets. Open hinge cover and remove pins. Set flexion first by moving hinge arms until an opening is seen through the desired flexion setting. Drop pins into desired setting, snapping them into place. Set both sides of the hinge (left/right) at the same degree of flexion. Set extension by moving hinge arms until an opening is seen through the desired extension setting. Drop pins into desired setting, snapping them into place. Set both sides of the hinge (left/right) to the same degree of extension. Close hinge cover. Insert hinges back into pockets.

### Caution! Please read carefully prior to product use:

#### CAUTION:

Wear as recommended. Always consult a physician when encountering persistent pain, discomfort, inflammation, or swelling.

#### WARNING:

Contact neoprene material with skin may cause irritation. Persons with a known susceptibility to dermatitis or with broken skin should avoid contact with neoprene. Anyone who develops skin irritation from contact with this material should discontinue its use immediately. Should not be worn by persons with phlebitis or circulatory disorders. Not recommended to be worn while sleeping.

#### CLEANING INSTRUCTIONS:

Fasten closures. Hand wash in cold water with mild soap. Wrap in a towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.



### Order information

Part Number: EZW33 Suggested HCPC: L1832/L1833 Size: Universal

#### How to order

Email: [orders@manamed.net](mailto:orders@manamed.net), Fax: 888-508-0712, Call: 888-508-0712

5240 W Charleston Blvd.  
Las Vegas NV 89146  
[www.ManaMed.com](http://www.ManaMed.com)



**MANAMED™**