

Royal Boot Air™

Pneumatic Walker Boot

- Rocker bottom to stimulate natural gait.
- Athletic bottom design.
- Large air bladder creates a larger surface for contact pressure to control swelling and ensure a snug fit.
 - To inflate, close round valve by turning it clockwise and add air pressure by pressing air bulb until desired fit.
 - To deflate, open round valve by turning counter-clockwise until air is released.
- Available sizes: Short or Tall; Extra Small through Extra-Extra Large.



Part No. SRB001**, RB001**, Suggested HCPC L4360, L4361

ManaEZ™ Ankle Stirrup

Ankle Support Brace

- Provides support and compression for ankle sprains and helps reduce swelling.
- Ideal for: acute injury, post-operation, ankle sprains, instability.
- Low profile for minimal bulk.
- Universal sizing.
- Adjustable straps allow for left or right usage.



Part No. EZAS4350, Suggested HCPC L4350

Ankle Kit

- ManaEZ™ Boot Air
- ManaFlexx 2™
- PlasmaFlow™
- ManaEZ™ Ankle Ice

Part No. AK001,
Suggested HCPC L1902,
L4361, E0745, E0676



ManaEZ™ Ankle Ice

- Quick ankle wrap orthosis with a removable ice pack.
- Inflation system that helps reduce pain and swelling.
- Universal size.



Part No. EZAnkle01, Suggested HCPC L1902



ManaFlow™

Pneumatic Compression Device

- Prescribed for:
 - Lymphedema
 - Chronic venous insufficiency
 - Prevention of DVT (Reduce edema and enhance venous return)
- Patient compliance and pressure monitor
- Four chambers; tubeless
- Portable and lightweight

Part No. MFLOW52 (Calibrated), Suggested HCPC E0652

Part No. MFLOW51 (Non-Calibrated), Suggested HCPC E0651

ManaEZ™ Boot Air

Pneumatic Walker Boot

- Rocker bottom to stimulate natural gait.
- Pneumatic pump;
 - To put air inside boot, turn valve clockwise until secure, then squeeze bulb until adequate pressure is felt.
 - To remove air from boot, turn valve counterclockwise.
- Large air bladder that creates a larger surface for contact pressure to help control swelling and ensure a snug fit.
- Available sizes: Short or Tall; Extra Small through Extra Large.

Part No. EZBA01**, Suggested HCPC L4360



PlasmaFlow™

Vascular Therapy Device

- Inflatable tube massage system that helps prevent Deep Vein Thrombosis by stimulating blood flow in the extremities (stimulating muscle contractions).
- Electronically-controlled pump delivers set amount of air to the leg cuffs, which then compress the calf or calves to aid in blood flow.
- Reduce post-operative pain and swelling.
- Reduce wound healing time.
- Portable; rechargeable battery.

Part No. PF0001, Suggested HCPC E0676



ManaFlexx 2™

Portable Electric Muscle Stimulator

- TENS and NMES;
 - Relaxes muscle spasms and prevents muscle atrophy.
 - Increases local blood circulation and healing.
 - Immediate post-surgical stimulation to prevent venous thrombosis.
 - Maintaining/increasing range of motion.
 - Temporarily relief of pain.
 - Symptomatic relief and management of chronic pain.

Part No. MF002-RX, Suggested HCPC E0745



ManaEZ™ Spring Crutch

Aluminum Crutch

- Universal sizing; Adjustable height.
- Comfortable and Lightweight.
- Absorbing spring at the bottom.
- Limited Warranty.

Part No. EZCRUTCH01SP,
Suggested HCPC E0117 x2



ManaMed believes in *Healing Through Innovation* and is dedicated to helping physicians treat and elevate their patient care by providing patients with high quality and innovative foot and ankle solutions to help them heal quicker.



To learn more about ManaMed, scan the QR code below



Manufactured By:



MANAMED™

5240 W Charleston Blvd.,
Las Vegas, NV 89146

888-508-0712 (toll free)
CustomerService@manamed.net
www.manamed.com

ManaMed and the ManaMed logo are trademarks of ManaMed, Inc.

Foot and Ankle Solutions



MANAMED™