

INSTRUCTIONS FOR USE

Mana **EZSLING**
Universal



INSTRUCTIONS FOR USE

FITTING INSTRUCTIONS

1. Place injured arm in sling.
2. Place thumb in thumb loop.
3. Place strap around neck and secure with hook and loop.
4. Place waist strap around waist and secure with hook and loop.

WASHING INSTRUCTIONS

1. Attach hook and loop straps together.
2. Hand or machine wash in mild detergent using cold water.
3. Air dry.

SUPPORTS + PROTECTS

Universal Size

Breathable, cotton material

Padded shoulder strap

Easily adjustable

Hook and Loop closure

Latex Free



MANAMED

ManaMed Inc.
5240 W Charleston Blvd.
Las Vegas NV 89146
United States
Ph: 888-508-0712
www.ManaMed.com

©2020 ManaMed Inc. All rights reserved.
ManaMed, the ManaMed logo,
are trademarks of ManaMed Inc.
Registered in the U.S.